

Development of a Next Generation Sequencing Panel to Assess Hereditary Cancer Risk that Includes Clinical Diagnostic Analysis of the *BRCA1* and *BRCA2* Genes

Benjamin Roa,¹ Karla Bowles,¹ Satish Bhatnagar,² Natalia Gutin,¹ Adam Murray,¹ Bryan Wardell,² Mark Bastian,² Jeffrey Mitchell,² Jian Chen,² Thanh Tran,² Deborah Williams,² Jenny Potter,² Srikanth Jammulapati,² Michael Perry,² Brian Morris,² and Kirsten Timms²
 1 - Myriad Genetic Laboratories, Inc., Salt Lake City, Utah 2 - Myriad Genetics, Inc., Salt Lake City, Utah


Introduction

- Approximately 7% of breast and 11-15% of ovarian cancers are estimated to be due to germline DNA mutations.
 - The majority of these mutations occur in the *BRCA1* and *BRCA2* genes, but mutations in additional genes have been shown to convey significant cancer risk.
- Sanger DNA sequencing has been the gold standard for detecting genetic variants.
 - Labor-intensive
 - Costly for the analysis of large gene panels
- Next Generation Sequencing (NGS) platforms could be another sequencing alternative.
 - Efficient analysis of larger gene panels
 - Lack of standardization of sample preparation, NGS platforms, and data analysis
- Optimized assay design and validation are critical to maximize the analytical sensitivity and specificity of NGS assays and to ensure high quality interpretation for clinical decision making.

Methods

- Initial assessment of analytical sensitivity and specificity was performed by comparing *BRCA1* and *BRCA2*.
- Testing performed on 1,864 anonymized patient samples which had previously undergone Sanger sequencing of *BRCA1* and *BRCA2*
- RainDance emulsion PCR system used to enrich targets for next generation sequencing
- Barcoded samples for multiple patients were pooled and loaded onto Illumina HiSeq and MiSeq sequencers
- NGS data analysis performed using a combination of commercial and laboratory-developed informatic tools

Results



Application - Development of 25 Gene Panel Assay

- Optimization of *BRCA1/BRCA2* analysis by NGS facilitated assay validation of 25-gene hereditary cancer panel that includes sequencing and large rearrangement analysis.
- BRCA1/BRCA2* testing is part of the 25 Gene NGS Panel
- We developed and validated a 25-gene NGS hereditary cancer panel (myRisk Hereditary Cancer™) that evaluates a broad number of hereditary cancer syndromes to help define patients' cancer risk and management options with a focus on eight primary cancer sites.

Cancer Genes and Associated Syndromes Targeted by 25 Gene Panel

Syndrome Name	Breast	Ovarian	Colon	Endometrial	Melanoma	Pancreas	Gastric	Prostate	Other cancer/Clinical parameters
<i>BRCA1</i>	BR	OV				PA		PR	
<i>BRCA2</i>	BR	OV			ME	PA		PR	
<i>MLH1</i>		OV	CO	EN		PA	GA		OC
<i>MSH2</i>		OV	CO	EN		PA	GA		OC
<i>MSH6</i>		OV	CO	EN		PA	GA		OC
<i>PMS2</i>		OV	CO	EN		PA	GA		OC
<i>EPCAM</i>		OV	CO	EN		PA	GA		OC
<i>APC</i>			CO			PA	GA		OC
<i>MUTYH</i>			CO						OC
<i>CDKN2A</i>					ME	PA			
<i>PALB2</i>	BR					PA			
<i>STK11</i>	BR	OV	CO	EN		PA	GA		OC
<i>PTEN</i>	BR		CO	EN					OC
<i>TP53</i>	BR	OV	CO	EN	ME	PA	GA	PR	OC
<i>CDHI</i>	BR		CO				GA		OC
<i>BMPRIA</i>			CO			PA	GA		OC
<i>SMAD4</i>			CO			PA	GA		OC
<i>ATM</i>	BR					PA			OC
<i>BARD1</i>	BR								OC
<i>BRIP1</i>	BR	OV							
<i>CDK4</i>					ME	PA			
<i>CHEK2</i>	BR		CO					PR	
<i>NBN</i>	BR							PR	
<i>RAD51C</i>	BR	OV							
<i>RAD51D</i>		OV							

Technical Components and Validation of 25 Gene Panel

Next Gen Sequencing	Large Rearrangements	Review and Reporting
<ul style="list-style-type: none"> RainDance emPCR prep LR-PCR for <i>PMS2</i>, <i>CHEK2</i> Illumina HiSeq, MiSeq NGS Sanger confirmation 	<ul style="list-style-type: none"> Microarray CGH <ul style="list-style-type: none"> Custom array Confirm LRs MLPA for <i>PMS2</i>, <i>CHEK2</i> <ul style="list-style-type: none"> Gene vs pseudogene LR vs gene conversion 	<ul style="list-style-type: none"> Informatics, LIMS workflow Variant classification Unusual case pathway Reporting

- 100% concordance** on NGS/Sanger parallel sequencing
 - 100 anonymized DNAs with 3923 variants in hereditary cancer panel genes
- Estimate >99.92% analytical sensitivity, specificity for gene panel
- Large Rearrangements (LR) validated with:
 - Microarray for 23 genes (212 DNAs with 51 genomic positive controls)
 - MLPA for *PMS2* and *CHEK2* (110 DNAs with 5 genomic positive controls)
 - Genomic positive controls supplemented with synthetic controls

Conclusions

- Preliminary analysis of the *BRCA1* and *BRCA2* genes has facilitated assay optimization and led to a comprehensive validation of all 25 genes in the NGS panel.
- Validation studies show that a NGS gene panel designed to meet rigorous quality standards provides clinically actionable results equivalent to those obtained from Sanger DNA sequencing analysis.

Presented at ASHG - October 24, 2013