# 25-GENE PANEL TESTING AND INTEGRATED RISK MANAGEMENT TOOL IMPACTS MEDICAL MANAGEMENT IN HEREDITARY CANCER SYNDROME EVALUATION



Lucy R. Langer, MD,<sup>1</sup> Laurie Korst, RN, BSN,<sup>2</sup> Larry J. Geier, MD,<sup>3</sup> Jennifer Saam, MS, LCGC, PhD,<sup>2</sup> Kelsey Moyes, MStat,<sup>2</sup> Brent Evans, MS,<sup>2</sup> Lee S. Schwartzberg, MD<sup>4</sup>

1Compass Oncology, Portland, OR; 2Myriad Genetic Laboratories, Inc., Salt Lake City, UT; 3Kansas City Cancer Ctr, Kansas City, MO; 4The West Clinic, Memphis, TN

### The West Clinic Advancing Cancer Care. Together. New Partners with Methodist and UTHSC

#### BACKGROUND

- The identification of patients with hereditary cancer syndromes such as Hereditary Breast and Ovarian Cancer (HBOC) or Lynch syndrome (LS) leads to clinical management changes. Using next generation sequencing, more comprehensive gene panels with greater sensitivity have been developed.
- Integrating personal and family cancer history identified during the screening process with genetic test results can offer refined management recommendations.
- The aim of this study was to determine the clinical utility of hereditary cancer panel testing.
- This presentation is focused on the breast, colon, endometrial and ovarian cancer medical management decisions made by health care providers based on gene panel results.

#### METHODS

#### **SCREENING**

• Patients were identified using criteria for HBOC or LS and tested using a 25-gene hereditary cancer panel.

#### **TESTING**

- Panel genetic testing included next generation sequencing and rearrangement analysis of 25 genes with cancer risk data:
  - BRCA1, BRCA2 MLH1, MSH2, MSH6, PMS2, EPCAM
  - APC, MUTYH CDKN2A, CDK4
  - PALB2, CHEK2, SMAD4, BMPR1A, STK11, TP53, CDH1, PTEN, and ATM
  - NBN, BARD1, BRIP1, RAD51C, and RAD51D
- Affected patients evaluated in this study included those with a personal history of breast, ovarian, colon, endometrial, gastric, pancreatic, melanoma, and prostate cancers.

#### **EVALUATION: PRE- AND POST-TEST SURVEYS**

- Recommendations from testing incorporated the genetic test result and a
  personalized cancer risk and management tool (myRisk management tool
  (MMT)) based on a patient's personal and family history and professional
  guidelines.
- Health care providers were surveyed for their management advice to the patient for breast (BC), ovarian (OC), endometrial (EC) and colorectal (CRC) cancer before and after testing.
- Data based on matched pre- and post-test surveys for 1111 patients.
- Interventions were ranked from least aggressive (surveillance) to most aggressive (surgery).

Least
<b>Aggressive</b>
71991C331VC

Surveillanc

noprevention


Surgery

Most Aggressive

#### Table 1. Clinical Characteristics

Characteristic	N=1111
Age at Testing (Yrs)	
Mean	50.0
Range	16 - 88
Gender	
Female	1077 (96.9%)
Male	34 (3.1%)
Ancestry	
Western/Northern Europe	631 (56.8%)
Central/Eastern Europe	79 (7.1%)
African	59 (5.3%)
Latin American/Caribbean	48 (4.3%)
Ashkenazi	25 (2.3%)
Asian	19 (1.7%)
Native American	9 (0.8%)
Neareast/Mideast	3 (0.3%)
Multiple Ancestries Indicated	88 (7.9%)
None Specified	150 (13.5%)
Personal History of Cancer	
Affected	516 (46.4%)
Unaffected	595 (53.6%)
NCCN Guideline Adherence	
2013 HBOC Only	734 (66.1%)
2012 Lynch syndrome Only	82 (7.4%)
Both HBOC and LS	214 (19.3%)
Neither	81 (7.3%)

#### Figure 1. Mutations Detected in Positive Patients




#### Table 2 Case Studies

 lable 2. Co	lable 2. Case Studies		
	Mutations Detecte		
Positive Pa	tients		
Patient #1	BRCA2		
Patient #2	CHEK2		
Patient #3	PMS2		
Patient #4	PMS2		
<b>Negative P</b>	atients		
Patient #5			
Patient #6			
Patient #7			

Patient #8

#### Figure 2. Information Used to Make Medical Management Decisions Table 3. Changes in Management Recommendations

\_ \_ \_



#### Bilateral mastectomy Breast MRI

	Positive Test Result (n=81)	Negative Test Result (n=1030)
More Aggressive	27 (33.3%)	78 (7.6%)
Mixed	13 (16.0%)	89 (8.6%)
Less Aggressive	5 (6.2%)	38 (3.7%)
Minor Change	9 (11.1%)	32 (3.1%)
No Change	18 (22.2%)	769 (74.7%)

Post-Test Recommendation

CA-125

Bilateral mastectomy

Pelvic exam

Pelvic exam

**Breast MRI** 

Transvaginal ultrasound

Colonoscopy

Transvaginal ultrasound

Clinical breast exam

Mammography

Colonoscopy

BSO

**Breast MRI** 

Pelvic exam

Pelvic exam

Colonoscopy

**Breast MRI** 

Colonoscopy

Transvaginal ultrasound

Chemoprevention

Mammography

Colonoscopy

Less

More

Less

More

Less

More

More

More

Less

More

None

More

None

More

More

More

More

Less

Less

Less

More

Less

\*Surveys for an additional nine (11.1%) positive test result patients and 24 (2.3%) negative test result patients indicated a change in medical management on the post-test survey but not the type of change.

#### RESULTS

BSO

Lumpectomy and Radiation

Transvaginal ultrasound

None

Clinical breast exam

Chemoprevention

None

Pelvic exam

None

Bilateral mastectomy

None

BSO

Mammography

Pelvic exam

None

None

Mammography

None

Chemoprevention

BSO

Bilateral mastectomy

None

Survey Category

OC

BC

CRC

EC

BC

CRC

OC

BC

OC

CRC

BC

CRC

CRC

BC

Yes

No

No

No

No

## Transvaginal US, CA-125, BSO, Ovarian chemoprevention, Mammography, Breast MRI, Breast chemoprevention, Clinical breast exam, Mastectomy, Pancreas cancer screening

Mammography, Breast MRI, Breast cancer reduction strategies, CRC screening

Colonoscopy, Endometrial sampling, Hysterectomy, BSO, Transvaginal US, CA-125, Upper endoscopy, Capsule endoscopy, Urinalysis, Pancreatic screening, Physical exam

Colonoscopy, Endometrial sampling, Hysterectomy, Mammography, BSO, Transvaginal US, CA-125, Upper endoscopy, Capsule endoscopy, Urinalysis, Pancreatic screening, Physical exam

Colonoscopy, Mammography, Breast MRI, Clinical breast exam

Mammography, Breast MRI, Clinical breast exam

Skin cancer screening

Colonoscopy, Mammography, Breast MRI, Clinical breast exam,
Breast cancer risk reduction strategies

#### CONCLUSIONS

- Over 74% of physicians used both the genetic test result and the personal and family cancer history from the MMT to make management decisions.
- Patients with positive genetic test results had management changes 77.8% of the time after testing.
- Patients with negative test results had management changes after testing in 25.3% of cases based on information gained from the MMT.