Prevalence of BRCA1/2-Negative Women Who Qualify for Adjunctive Breast MRI Screening

Heidi Gorringe, MS, CGC, Eric Rosenthal, PhD, Michael Van Orman, BS, John Kidd, MS, Krystal Brown, PhD, Royce T. Adkins, MD

Myriad Genetic Laboratories, Inc., Salt Lake City, UT

BACKGROUND

- The National Comprehensive Cancer Network (NCCN) and the American Cancer Society (ACS) recommend breast MRI in addition to mammography for women with an estimated lifetime risk of breast cancer >20%.1,2
- Originally, this recommendation applied to women identified using breast cancer risk models "largely dependent on family history", and women identified with pathogenic variants (PVs) in the high-penetrance breast cancer risk genes: BRCA1, BRCA2, PTEN, TP53, CDH1 and STK11.
- Recently, NCCN expanded the breast MRI recommendation to women identified as carrying PVs in the genes ATM, CHEK2 and PALB2,1 based on the estimated >20% lifetime breast cancer risk associated with these genes.

OBJECTIVE

- This analysis compares the proportion of women identified as having a >20% lifetime risk for breast cancer as an outcome of genetic testing for BRCA1 and BRCA2 versus expanded testing utilizing a panel that includes the genes ATM, CHEK2 and PALB2.
- We also determined the percentage of women identified with PVs in ATM, CHEK2 and PALB2 who would have been identified as having a >20% risk based on family history assessement

Table 2. Breakdown of ATM, CHEK2 and PALB2 Pathogenic Variants and Claus Assessment for Female Carriers

Gene	Observations	Patients with Breast Cancer	Patients with estimated risk >20% based on family history*
ATM	315	167 (53.0%)	63 (20.0%)
CHEK2	379	185 (48.8%)	80 (21.1%)
PALB2	229	147 (64.2%)	63 (27.5%)
ATM & CHEK2	2	0	2
ATM & PALB2	2	0	1
CHEK2 & PALB2	5	5	0
Total	932	504 (54.1%)	209 (22.4%)

Note: Percentages were calculated according to gene. Percentages were not calculated for rows with fewer than 10 patients.

* Based on Claus model, without including personal diagnosis of tested patient

METHODS

- All data were derived from clinical testing ordered for 45,661 female patients over a 1 year period from September 2013 to November 2014.
- Women were tested using a 25-gene hereditary cancer panel that includes BRCA1, BRCA2, ATM, CHEK2, and PALB2 (Table 1).
- Testing included full sequencing and comprehensive large rearrangement analysis, except for EPCAM (large rearrangement only).
- PVs are those that received a laboratory classification of Deleterious or Likely Deleterious.

Table 1. 25-Gene Hereditary Cancer Panel

Genes	Associated Cancers		
BRCA1, BRCA2	Breast, Ovarian, Pancreatic, Prostate		
ATM	Breast, Pancreatic		
CHEK2	Breast		
PALB2	Breast, Pancreatic		

Other Panel Genes


Genetic Testing

APC, BARD1, BMPR1A, BRIP1, CDH1, CDK4, CDKN2A, EPCAM, MLH1, MSH2, MSH6, MUTYH (biallelic), NBN, PMS2, PTEN, RAD51C, RAD51D, SMAD4, STK11, TP53

Family History Assessment

- All clinical information was obtained from the test request form.
- The Claus model was applied to women carrying PVs in ATM, CHEK2, and PALB2 to determine which patients would have been identified as having a >20% risk of breast cancer based on family history alone.
- The Claus tables provide a lifetime risk, defined as to age 80, based on the breast cancer history in a woman's first and second-degree relatives.3
- Affected individuals were assigned a risk estimate without consideration of existing diagnoses in order to determine if they would have been identified for more aggressive screening prior to their breast cancer.
- For the purposes of this analysis, the Claus tables were used to estimate a lifetime risk for a woman at age 29, without adjusting for her current age. This is analogous to how risk estimates are applied for patients with risk due to PVs in inherited breast cancer genes, which are not typically adjusted for
- This results in an over-estimate of women identified as having a >20% risk based on family history, as risk estimates based on the tables decline as women

Figure 1. Proportion of Women with a PV in ATM, CHEK2, or PALB2 Who Would Have Been Candidates for Breast MRI Based on Family History Assessment




Excludes women who also carried a PV in BRCA1/2; Based on Claus model, does not include personal diagnosis of patient

Presented at NCoBC - April 9-11, 2016

Met 20% Risk Threshold Based on Family History Alone

RESULTS

- 932 women were identified as carrying PVs in only ATM, CHEK2, or PALB2 (Table 2).
 - 9 women carried PVs in 2 of these genes.
 - 16 patients with a second PV in BRCA1 or BRCA2 were excluded from analysis.
- Only 209 (22.4%) of these women reached the >20% lifetime breast cancer risk threshold using the Claus tables (Figure 1).
 - This demonstrates that family history alone will not identify the majority of women with PVs in ATM, CHEK2, and PALB2.
- The percentage of unaffected female PV carriers meeting the >20% risk threshold (25.7%) was higher than that for affected patients (19.6%) (Figure 1).
 - This is consistent with the presumption that unaffected patients would require a stronger family history to be considered appropriate candidates for testing.
- 1,467 women were identified as carrying a PV in BRCA1 or BRCA2.
- The use of panel testing to detect PVs in ATM, CHEK2, and PALB2 resulted in a 63.5% increase over the total found in BRCA1 and BRCA2 (Figure 2).



CONCLUSIONS

- Relative to BRCA1 and BRCA2 testing alone, we found that clinical testing with a hereditary cancer risk panel that includes ATM, CHEK2, and PALB2 in this cohort resulted in a 63.5% increase in the identification of women who are eligible for MRI screening.
- Close to 80% of the women with mutations in ATM, CHEK2 and PALB2 would not have been identified as having a >20% lifetime breast cancer risk based on family history assessment alone.

REFERENCES

- 1. Daly M, Pilarski R, Axilbund JE, et al. Genetic/Familial High-Risk Assessment: Breast and Ovarian. Version 2.2016. NCCN Clinical Practice Guidelines in Oncology. 2016. http://www.nccn.org/professionals/physician_gls/pdf/genetics_screening.pdf.
- 2. American Cancer Society. Breast Cancer Prevention and Early Detection. 2015. http://www.cancer.org/acs/groups/cid/documents/webcontent/003165-pdf.pdf
- 3. Claus E. B., Risch N., Thompson W. D. *Cancer*, 1994; 73: 643-651.