

CCP SCORE IS A STRONG PREDICTOR OF OUTCOME IN SEVERAL PROSTATE CANCER COHORTS



MYRIAD

Steven Stone¹, Jack Cuzick², Dan Berney², Julia Reid¹, David Mesher², Gabrielle Fisher²,
Jerry Lanchbury¹, Alexander Gutin¹, Greg Swanson³ and on behalf of the Transatlantic Prostate Group

¹Myriad Genetics, Salt Lake City, UT; ²Cancer Research UK, London, UK; ³University of Texas Health Science Center, San Antonio, TX



BACKGROUND

- The natural history of newly diagnosed prostate cancer is highly variable and difficult to predict so improved tools are needed to more appropriately match treatment to a patient’s risk of progression.
- Previous data has shown that a 46-gene cell cycle progression (CCP) RNA signature is a robust predictor of biochemical recurrence after radical prostatectomy, and of prostate cancer-specific death in a conservatively managed cohort diagnosed by transurethral resection of the prostate.¹
- Here we report the prognostic utility of the CCP score obtained in 3 different clinical settings including needle biopsies in newly diagnosed men.

METHODS

- mRNA was extracted from formalin fixed and paraffin embedded (FFPE) tumor sections from 1) 366 U.S. patients after radical prostatectomy; 2) 337 conservatively managed (i.e. watchful waiting) UK patients diagnosed by TURP; and 3) 349 conservatively managed UK patients diagnosed by needle biopsy (Table 1).
- RNA levels of 31 CCP genes and 15 housekeeper genes were determined and a mean composite score calculated (CCP score).
- Clinical variables for multivariate analysis included Gleason score, baseline PSA, age, and stage.
- Primary endpoint was death from prostate cancer in the TURP and needle biopsy cohorts, and biochemical recurrence in the RP cohort.

RESULTS

- 1) CCP score was highly prognostic of outcome in all tested clinical settings (Table 2 and Table 3).
- 2) CCP score was highly predictive of biochemical recurrence after prostatectomy (p-value < 10⁻⁸; HR = 1.89). After adjustment for clinical parameters including PSA, Gleason, pStage and margins CCP score remained highly significant (p-value < 10⁻⁵; HR = 1.77).
- 3) In the TURP cohort, the CCP score was highly predictive of disease-specific death (p-value < 10⁻²¹; HR = 2.92). It remained a significant predictor of cancer death after adjustment for Gleason grade, PSA, Ki67 status, and cancer extent (p-value < 10⁻⁷; HR = 2.56).
- 4) In the needle cohort, CCP score was the strongest univariate predictor of cancer death (p-value < 10⁻⁹; HR = 2.02). It remained significant after adjustment for Gleason grade, PSA, Ki67, and extent of disease (p-value < 10⁻⁴; HR = 1.65).

CONCLUSIONS

- An mRNA expression signature based on CCP gene expression is prognostic in prostate cancer patients at diagnosis and after prostatectomy.
- CCP score provides important prognostic information that is not provided by other clinical or pathological variables.
- The CCP signature should be a valuable addition to clinical variables for differentiating aggressive from indolent disease.

Table 1. Patient characteristics for all 3 cohorts.
Numbers are median (IQR) or n (%) as appropriate.

Cohort	Post-RP ¹	TURP ¹	Needle biopsy ²
	N=366	N=337	N=349
Outcome	Time to biochemical recurrence	Death from prostate cancer	Death from prostate cancer
Events	138 (38)	75 (22)	90 (26)
Follow-up	9.5 (6.8, 11.0)	10.3 (5.9, 12.3)	10.3 (5.5, 11.6)
Age (years)	68 (63, 72)	71 (67, 73)	71 (66, 73)
Gleason Score			
<7	240 (66)	172 (51)	106 (30)
7	110 (30)	73 (22)	152 (44)
>7	16 (4)	92 (27)	91 (26)
PSA (ng/ml)	6.9 (4.5, 10.7)	8.3 (2.8, 21.0)	21.4 (11.9, 42.0)

Table 2. Summary of Cox proportional hazards univariate analysis.
Numbers are hazard ratio (95% CI) Hazard ratio for CCP score is for an increase in one score unit.

Cohort	Post-RP ¹		TURP ¹		Needle biopsy	
Outcome	Time to biochemical recurrence		Death from prostate cancer		Death from prostate cancer	
	Hazard Ratio (95% CI)	p-value	Hazard Ratio (95% CI)	p-value	Hazard Ratio (95% CI)	p-value
CCP score	1.89 (1.54, 2.31)	5.6 x 10 ⁻⁹	2.92 (2.38, 3.57)	6.1 x 10 ⁻²²	2.02 (1.62, 2.53)	8.6 x 10 ⁻¹⁰
Gleason score						
<7	1 (ref)	1.5 x 10 ⁻¹³	1 (ref)	3.7 x 10 ⁻¹⁹	1 (ref)	1.6 x 10 ⁻⁹
7	2.81 (2.01, 3.94)		5.20 (2.40, 11.29)		2.17 (1.16, 4.07)	
>7	6.32 (3.65, 10.93)		13.67 (6.90, 27.11)		5.85 (3.11, 11.01)	
log(1 + PSA) ng/ml	3.07 (2.38, 3.97)	3.4 x 10 ⁻¹⁷	2.30 (1.83, 2.88)	3.4 x 10 ⁻¹⁴	1.70 (1.31, 2.20)	4.2 x 10 ⁻⁵

Table 3. Summary of Cox proportional hazards multivariate analysis.
HRs are given per unit increase in CCP score.

Cohort	HR (95% CI)	p-value	Clinical variables included in model
Post-RP	1.74 (1.39 - 2.17)	3.3 x 10 ⁻⁶	Gleason, PSA, stage, margins
TURP	2.56 (1.85 - 3.53)	1.3 x 10 ⁻⁸	Gleason, PSA, Ki67, % positive chips
Needle biopsy	1.65 (1.31, 2.09)	3.0 x 10 ⁻⁵	Gleason, PSA, stage, hormone use, age

Figure 1. Risk of 10-year prostate cancer mortality by Combined risk (CCP + Gleason + PSA) for patients in the needle biopsy cohort.





Figure 2. 10-year predicted risk in needle cohort for model with CCP score compared to model with PSA and Gleason.
Blue = Gleason 6; Black = Gleason 7; and Turquoise = Gleason 8-10.



REFERENCES

1. Cuzick et al. *Lancet Oncol* 2011; 12(3): 245-55.