

# THE ROAD TO MANAGING RHEUMATOID ARTHRITIS

USE VECTRA® DA TO HELP GUIDE THE WAY

RA JOURNEY →

1.5M  
AMERICANS  
LIVE WITH RA

1

RECEIVE RA  
DIAGNOSIS

2

START  
MEDICINE

VECTRA DA  
SCORE  
**55**

Vectra DA is an **advanced blood test** that helps doctors and patients better understand RA disease activity.

**CAUTION**

Single biomarker tests may lack complete information<sup>1</sup>.

3

ASSESS  
RESPONSE

12 biomarkers are combined into a **single Vectra DA score** indicating the current level of disease activity

VECTRA DA  
SCORE

LEVEL OF  
DISEASE  
ACTIVITY

45 to 100  
30 to 44  
1 to 29

HIGH  
MODERATE  
LOW

4

ADJUST  
MEDICINE

5

MONITOR  
WELL-CONTROLLED RA

VECTRA DA helps to assess:

- Disease activity
- Risk of joint damage
- Response to treatment

VECTRA DA  
SCORE  
**45**

VECTRA DA  
SCORE  
**18**

## CHANGE THE WAY YOU LOOK AT RA WITH A SINGLE SCORE

Vectra DA measures **12 key proteins** that represent **multiple biological pathways** linked to RA, providing:

- A more complete picture of RA disease activity than single biomarkers such as erythrocyte sedimentation rate (ESR) and C-reactive protein (CRP)<sup>2</sup>
- A better assessment of risk of joint damage than either ESR or CRP<sup>2,3</sup>

With a **single score** that can easily be tracked over time, Vectra DA offers doctors and patients a **more comprehensive understanding** of RA disease activity that can be used to make **more informed treatment decisions**.